



1. Rapid and Deep Spatial Phenotyping

Single-cell spatial phenotyping has transformed cancer research and is poised to play formative roles in the development of effective therapeutic strategies. Here, we present ultrahigh-plex single-cell spatial phenotyping of whole-slide human FFPE tissues with 100+ protein biomarkers encompassing *immune cell lineage, activation states, immune checkpoints, tissue structure, apoptosis, DNA damage & metabolism.*



Deployment of our 103-plex deep spatial phenotyping panel occurs with rapid turn-around time on the PhenoCycler™-Fusion, an automated spatial biology platform for ultrahigh-plex imaging. The PhenoCycler (formerly CODEX®) integrates seamlessly with the Fusion for end-to-end spatial imaging at high resolution and scale.


2. Design and Development of 103-plex Panel




The design of the 103-plex antibody panel is based on the hallmarks of cancer. The panel includes markers for cell lineage, immune activation and checkpoints, cellular energetics, mediators of proliferation and metastasis and stress responses, and more. Each marker and each module has been carefully selected to reveal unique information on different pathways and, when multiplexed together, provides an integrated overview of the landscape of cancer progression

3. Largest Single-cell Spatial Phenotyping Dataset from a Single Sample

3.1 Deploying a 103-plex Antibody Panel in FFPE Tissue




3.2 Cross-Correlational Analysis of 103 Markers Across the Entire Sample



3.1 Whole-Slide Imaging at Single-Cell Resolution (left; biomarkers as indicated) and Spatial Phenotyping (right) of a human FFPE tonsil into 11 distinct cell types. 3.2 Cross-correlation Analysis of 103 markers indicates cell lineages and functional niches with varying metabolic, immune, and stress signatures.


3.3 Single-Cell Phenotypic Characterization of 2.25 Million Cells



3.3 Single-Cell Phenotyping of 2'253'539 cells in a human tonsil shows an abundance of B- and T-cells. The heatmap (with markers listed on the right) shows a curated clustering dendrogram with cell types (columns) corresponding to the legend as indicated. The UMAP plot (top left) shows distinct phenotype clusters sorted by color and geography. The abundance of each cell type is quantified in the accompanying pie chart (bottom left).


4. Ultra High-Plex Single-Cell Spatial Phenotyping of Human FFPE Oropharyngeal Squamous Cell Carcinoma

4.1 Single-Cell Spatial Phenotyping of Human Oropharyngeal Squamous Cell Carcinoma reveals 14 Distinct Cell Types



4.1 Whole-Slide Imaging at Single-Cell Resolution of a FFPE human head and neck tumor (left; biomarkers as indicated) and Spatial Phenotyping (right) of all 835'823 cells into 14 distinct cell types. The heatmap (far right) includes a curated clustering dendrogram with cell types corresponding to the legend as indicated below it. The UMAP plot and the pie chart show the abundance of distinct phenotype clusters sorted by color and geography. 4.2 Deep Spatial Phenotyping shows the complex topography and disruption of normal tonsil architecture in tumor tissue (scatter plot; far left). Further analysis reveals 4 distinct tumors within the tissue (shown here in blue, purple, green and orange) with varying abundance of immune, proliferating, epithelial, and vascular cell types compared to the normal lymphoid area (red) and the normal squamous epithelium (black). The accompanying bar chart demonstrates the intra-tumor heterogeneity and highlights the importance of single-cell spatial phenotyping to understand the complexities and nuances of tumor progression.

4.3 The 103-Plex "Cancer Hallmark" Panel reveals Distinct Immune, Metabolic, and Stress Signatures reflecting Intra-tumor Heterogeneity



4.3 Metabolic, Stress and Immune Signatures from Ultra High-Plex Spatial Phenotyping. The heatmaps (left) show a curated clustering dendrogram of metabolic, apoptotic, stress, and immune markers within CD4+ T cells, CD8+ T cells, M1 macrophages, M2 macrophages, and tumor cells. Intra-Tumor Heterogeneity is shown by multiplex immunofluorescence images (top) of CD68 (red), a macrophage marker, and CD8 (green), a cytotoxic T cell marker demonstrating the differences in immune infiltration between the different tumors. Differential expression of Glucose-6-phosphate dehydrogenase (red), a key enzyme in glucose metabolism, and BAK (green), a pro-apoptotic regulator reflects the unique metabolic and stress features within the tumor.

5. The Need for Deep, Unbiased Spatial Analysis at Scale

This study demonstrates the power of rapid, deep single-cell spatial phenotyping enabled by the PhenoCycler-Fusion system. We present a meticulously designed, ultra high-plex antibody panel to decipher critical mechanistic insights into tumor biology. Whole-slide, single-cell spatial phenotyping uncovers the intricacies of intra-tumor heterogeneity and can aid in revealing mechanisms underlying clinical response and therapeutic resistance.