

## 1 Background

We describe two advances in multispectral fluorescence immunohistochemistry (fIHC), a powerful tool for quantifying interactions within the tumor microenvironment.


1. A **fully-automated 8-plex, 9-color assay plus DAPI counterstain** on the same tissue section.
2. A novel scanning method that produces a **multispectral whole slide scan of 6 markers plus DAPI counterstain in ~6 minutes** (1x1.5 cm tissue section).

## 2 Methods

FFPE samples of primary tumors were immunostained using Opal™ reagents manually or on a Leica BOND RX™. Imagery was acquired on a Vectra Polaris® automated imaging system and analyzed with inForm®, MATLAB®, and R software.

### Multiplex Staining with Opal Reagents


Opal™ reagents allow multiplex fIHC staining with signal amplification and any combination of mouse and/or rabbit primary antibodies.



**Fig 1. Opal™ Detection.** The Opal Polymer HRP amplifies IHC detection by covalently depositing multiple Opal fluorophores near the detected antigen. Then, antibodies are stripped to allow for sequential labeling of multiple markers.


### Multispectral Imaging on Vectra Polaris

**Fig. 2. Multispectral imaging on the Vectra Polaris is built upon an epifluorescence light path (below, left).** Different combinations of agile LED bands, bandpass excitation filters, bandpass emission filters, and a liquid crystal tunable filter (LCTF) are used to select narrow spectral bands that reach the imaging sensor.




For each spectral band, an image is acquired and added to a 'data cube' that contains up to 40 spectral layers (above, right). The data from all spectral layers is then linearly unmixed using previously-determined pure emission spectra for each fluorophore using inForm® software. Intensity values in the resulting 'unmixed' image are directly related to the amount of each dye present.


## 3 Results: 9-color Multispectral Imaging



**Fig. 3. 8-plex, 9-color panel on human lung cancer section.** Two new Opal™ reagents (Opal Polaris 480 and Opal Polaris 780) were combined with the currently available Opal 7-Color Automation IHC Kit to stain and distinguish 8 markers plus DAPI when imaged on the Vectra Polaris®: CD20 (Opal Polaris 480), PD-L1 (Opal 520), CD8 (Opal 540), FoxP3 (Opal 570), CD68 (Opal 620), PD-1 (Opal 650), Ki67 (Opal 690), and PanCK (Opal Polaris 780). Colors assigned to each marker, and associated component planes, are shown on the right.



**Fig 4. Phenotype densities from 9-color panel.** **A**) Whole slide view of lung cancer section shown above in Fig. 3. Multispectral fields were acquired at the locations indicated by the white boxes. **B**) Average densities of each phenotype marker within the acquired fields are shown for each section. Tissue was categorized into tumor and stroma regions using CK positivity by a trained segmentation algorithm in inForm® software. PD-L1 positivity was determined by a fixed intensity threshold; all other phenotypes were determined by a trained inForm® phenotyping algorithm.



**Fig 5. Cell proliferation (Ki67+) and PD-L1/PD-1 proximity assessment.** The 8 immune markers combine to generate more than 20 phenotypes relevant to immunology. A subset of these markers were studied in relation to proliferation state (**A**) and local PD-L1 proximity (**B**) to characterize the tumor microenvironment. For example, (**A**) shows that CD8+ and CD68+ cells were more likely to be proliferating (Ki67+) in tumor vs. stroma, with differences in proliferation observed in the presence or absence of PD-1 or PD-L1. (**B**) shows the high proximity between PD-1+ and PD-L1+ cells in this sample.

## 4 Results: 7-color Multispectral Whole Slide Scans


Field-based multispectral imaging workflows can accommodate a wide range of fluorophores and up to 9 colors, but can be time consuming as they require up to 50 spectral layers to unmix 9 fluorophores, and often require exposure times in the hundreds of milliseconds.

We have developed complementary high-throughput multispectral scanning approach by optimizing a multispectral workflow for a specific set of 7 fluorophores.

High-throughput multispectral scanning and unmixing performed comparably to field-based multispectral imaging, and outperformed conventional scanning by:

- **Reducing autofluorescence contributions** for all immune markers, **lowering the limit of detection and extending the dynamic range** of some channels by more than 30-fold.
- **Reducing crosstalk** from more than 8% to under 3% (typically <0.5%), thereby **increasing signal accuracy** and reducing false colocalization between non-colocalized markers.

For more details on 7-color whole slide MSI, see related posters via QR code at top right.



**Fig 6. Cell density and interaction density across the whole slide.** **A**) Whole slide MSI of human lung cancer section captured in 6 minutes, shown as composite image with marker colors indicated in key. Cells were phenotyped in inForm®, and interactions assessed with R and phenoptr. **B**) Density contours of CK+ (left), CD8+ (middle), and CK+ within 30 μm of a CD8+ cell (right). **Bottom**) Zoomed in views of A illustrate differences in CD8+ T-cell (yellow) infiltration within the tissue.

## 5 Conclusions

We introduce a **9-color fIHC assay** that distinguishes 8 markers plus DAPI counterstain on the same tissue section, increasing the depth of cellular interactions that can be studied within the tumor microenvironment.

Additionally, we introduce a **whole slide multispectral imaging** method that provides rich quantitation of interactions among 6 markers at length scales spanning from cell biology to tumor physiology.